

CD4 and CD8 T cell-specific immune responsiveness of QFT[®]-Plus

Background

TB-specific CD8+ T cell (CD8) release of IFN-g are believed to a biomarker of intracellular TB bacilli burden. This assumption is based on immunology studies that have found TB-specific CD8 responses more frequently in those with active TB disease compared to latent infection (1, 2), associated with recent exposure to TB (3), in active TB subjects with HIV co-infection and who are young children (4, 5), and declining responses in patients on anti-tuberculosis treatment (6).

Unlike QFT-TB Gold In Tube (QFT-GIT), QuantiFERON[®] TB Gold Plus (QFT[®]-Plus) is designed to utilize both CD4+ T cell (CD4) and CD8 specific immune responses. Additionally, the isolated CD8 response can be quantified by subtracting the response of QFT-Plus' TB1 tube (TB1) from TB2 tube (TB2) to provide immunologic CD8 information that may aid in the patient's clinical assessment and management. In this report, I will introduce the recent publication on validating the CD4 and CD8 immune responses of the TB1 and TB2 tubes of QFT-Plus.


New Publication

Allen et al. studied QFT-Plus TB1 and TB2 immune responses by using CD4 and CD8 clones (7).

A CD4 clone (A. 38.1 # 1, recognizing the long-chain peptide of CFP-10) and three CD8 clones (B. D481 C10, C. D432 H12, D. D480 F6, CFP-10 recognition of short peptides) were used for their study. Each of the clones was re-suspended in 10% human serum RPMI medium with various concentrations (10^3 / mL = K / mL), and 1 mL of cell suspension was added to both TB 1 and TB 2 tubes, then mixed and incubated at 37 °C. After incubation, each tube was centrifuged, the IFN-g value (IU/mL) in supernatant was measured by QFT-Plus ELISA. The results are shown in Figure 1 and Table 1.

Immune responses of CD4 clone (A) were similar in TB1 and TB2 in dose-response manner, resulting in increased immune responsiveness (Figure 1A). The immune response of one CD8 clone (B) to TB1 was slightly elevated in response to the number of clones but the responsiveness to TB2 was overwhelmingly higher than that to TB1 (Figure 1B). Furthermore, responsiveness of CD8 clones (C) and (D) to TB1, was slight, and responsiveness to TB2 was more than 100 times responsive compared to TB1 (Table 1-C, D).

The results of this study that utilized specific CD4 and CD8 clones validate that the antigens of TB1 is CD4 specific and TB2 is specific to both CD4 and CD8, as QFT-Plus was designed.



CD4 clone
D16038.1#1 No. of cells K/mL

CD8 clone
D481 C10 No. of cells K/mL

Figure 1. Immune responses of CD4 and CD8 clones to QFT-Plus TB1 and TB2 (A and B)

Table 1. Immune responses of CD8 clones to QFT-Plus TB1 and TB2 (C and D)

	QFT-Plus Blood Collection Tubes	No. of cells 20K	IFN- γ IU/mL			
			40K	60K	80K	100K
C.						
CD8 clone	TB1	0.12	0.24	0.22	0.33	0.51
Cd432H12	TB2	36	73	149	254	266
D.						
CD8 clone	TB1	N/D	N/D	N/D	0.53	0.98
D480F6	TB2	N/D	N/D	N/D	179	149

Note: If more than 10 IU/ml, the measured value was calculated by diluting in 10 to 100 times.

N/D: not done, K=10³ No. of cells.

Comment: This study provides direct evidence of the specificity of CD8 T-cell mediated response in QFT-Plus TB2 tubes. Compared with conventional existing IGRA in immunocompromised cases, QFT-Plus has been reported to be more sensitive compared to QFT-GIT due to CD8 signal in addition to CD4 (8–11). Further, the TB2-TB1 value is represented as the isolated immune response of CD8. The TB2-TB1 value of QFT-Plus suggests the possibility of additional clinical

significance, such as a CD8 immune response associated with a higher intracellular TB condition such as active TB, developing active tuberculosis (12), and recent infection in latent tuberculosis infection (13). The potential clinical significance of the CD8 response of QFT-Plus shall be clarified with the accumulation of data obtained through future and ongoing studies.

Mr. QFT

References

1. Day, C.L. et al. (2011) Functional Capacity of Mycobacterium tuberculosis-Specific T Cell Responses in Humans Is Associated with Mycobacterial Load. *J. Immunol.* **187**, 2222-2232
2. Rozot, V. et al. (2013) Mycobacterium tuberculosis-specific CD8+ T cells are functionally and phenotypically different between latent infection and active disease. *Eur J Immunol.* 2013 Jun; **43**(6):1568-77
3. Nikolova, M. et al. (2013) Antigen-specific CD4- and CD8-positive signatures in different phases of Mycobacterium tuberculosis infection. *Diagn. Microbiol. Infect. Dis.* **75**, 277-281.
4. Chiacchio, T. et al. (2014) Polyfunctional T-cells and effector memory phenotype are associated with active TB in HIV-infected patients. *J. Infect. Dis.* **69**(6):533-545
5. Lanicioni, C. et al. (2012) CD8+ T cells provide an immunologic signature of tuberculosis in young children. *Am J Respir Crit Care Med.* Jan 15; **185**(2):206-212.
6. Nyendak M. Et al. (2014) Mycobacterium tuberculosis Specific CD8+ T Cells Rapidly Decline with Antituberculosis Treatment. *PLoS ONE* **8**, e81564. Epub.
7. Allen et al. (2018)– Tuberculosis (2018), Characterization of specific CD4 and CD8 T-cell responses in QuantiFERON TB Gold-Plus TB1 and TB2 tubes, *Tuberculosis* **113**, 239–241
8. Telisinghe et al. (2017) The sensitivity of the QuantiFERON®-TB Gold Plus assay in Zambian adults with active tuberculosis. *Int. J. Tuberc. Lung Dis.* **21**(6), 690–696.
9. Igari, H, et al. (2018) Lymphocyte subset analysis in QuantiFERON-TB Gold Plus and T-Spot.TB for latent tuberculosis infection in rheumatoid arthritis. *J Infect Chemother.* Feb. **24**(2), 110-116.
10. Chien, J. et al. (2018) QuantiFeron-TB Gold Plus is a More Sensitive Screening Tool than QuantiFeron-TB Gold In-Tube for Latent Tuberculosis Infection among Older Adults in Long-Term Care Facility. *J. Clin. Microbiol.* Doi. 10.1128/JCM.00427-18
11. Fukushima, K at al. (2018) Comparison Study of Sensitivity of QuantiFERON TB Gold Plus with Existing IGRAs in the Patients with Active Pulmonary Tuberculosis, *Kekkaku*, Vol.93 No.10, 517-523
12. Barcellini L et al. (2016) First evaluation of QuantiFERON-TB Gold Plus performance in contact screening. *Eur Respir J.* Nov **48**(5):1411-1419
13. Pieterman ED et al. (2018) A multicentre verification study of the QuantiFERON®-TB Gold Plus assay. *Tuberculosis* **108** 136–142

QFT-Plus is an in vitro diagnostic aid for detection of Mycobacterium tuberculosis infection (including disease) and is intended for use in conjunction with risk assessment, radiography, and other medical and diagnostic evaluations. QFT-Plus results alone cannot distinguish active TB disease from latent infection. QFT-Plus Package Inserts, available in multiple languages, as well as up-to-date licensing information and product-specific disclaimers can be found at www.QuantiFERON.com.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Sample to Insight®, QFT®, QuantiFERON® (QIAGEN Group) 2400909 02/2019 © 2019 QIAGEN, all rights reserved.

インターフェロニンγ遊離試験キット
QuantiferON TB ゴールド プラス

保険適用 体外診断用医薬品
製造販売承認番号：23000EZ00004000

【製造販売業者】
株式会社 キアゲン

【お問い合わせ先】
株式会社 キアゲン カスタマーサポート

〒104-0054 | 東京都中央区勝どき3-13-1 | Forefront Tower II
Tel:03-6890-7300 | Fax:03-5547-0818
www.QuantiFERON.com

真空密封型採血管
QuantiferON TB ゴールド プラス チューブ

管理医療機器
認証番号：229AFBZX0004000